



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
Telephone: 901-382-8716  
Fax: 901-333-8223  
Email: info@meridianlifescience.com  
www.MeridianLifeScience.com

### SPECIFICATION SHEET

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** C86542M **Lot #:** 7J29707

**Description:** MAb to B. anthracis (LF)  
Monoclonal Antibody to Bacillus anthracis Lethal Factor (Anthrax)

**Specificity:** Bacillus anthracis Lethal Factor

**Clone:** B743M

**Host Animal:** Mouse. Hybridization of Sp2/0 myeloma cells with spleen cells from BALB/c mice. **Isotype:** IgG<sub>2a</sub>

**Source:** Ascites

**Immunogen:** Highly purified Bacillus anthracis lethal factor

**Format:** Purified, Liquid

**Purification:** ≥95% pure (SDS-PAGE). Protein G chromatography

**Concentration:** 1.4mg/ml (OD280nm, E<sup>0.1%</sup>=1.4)

**Affinity Constant:** Not determined

**Buffer:** PBS, pH 7.4

**Preservative:** 0.1% Sodium azide

**Applications:** Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended pair for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
C86543M	C86542M

**Storage:** Store at 2–8°C

**Warning:** This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**