



# Meridian

## Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road  
Memphis, TN 38134  
Telephone: 901-382-8716  
Fax: 901-333-8223  
Email: info@meridianlifescience.com  
www.MeridianLifeScience.com

### SPECIFICATION SHEET

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

**Catalog #:** H01289M **Lot #:** n/a

**Description:** MAb to Adiponectin  
Monoclonal Antibody to Human Adiponectin

**Specificity:** Human serum Adiponectin

**Clone:** B863M

**Host Animal:** Mouse. Hybridization of Sp2/0 myeloma cells with spleen cells from BALB/c mice. **Isotype:** IgG<sub>1</sub>

**Source:** Ascites

**Immunogen:** Native human adiponectin

**Format:** Purified, Liquid

**Purification:** >90% pure (SDS-PAGE). Protein A chromatography

**Concentration:** Lot specific

**Affinity Constant:** Not determined

**Buffer:** PBS, pH 7.4

**Preservative:** 0.1% Sodium azide

**Applications:** Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antibody pair for detection of low molecular weight Adiponectin in sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
H01289M	H01262M

**Storage:** Store at 2-8°C.

**Warning:** This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**