



Meridian

Life Science,® Inc.

Innovative Solutions. Trusted Partner.®

5171 Wilfong Road
Memphis, TN 38134
Telephone: 901-382-8716
Fax: 901-333-8223
Email: info@meridianlifescience.com
www.MeridianLifeScience.com

SPECIFICATION SHEET

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Catalog #: H01290M **Lot #:** n/a

Description: MAb to Adiponectin
Monoclonal Antibody to Human Adiponectin

Specificity: Human serum Adiponectin

Clone: B864M

Host Animal: Mouse. Hybridization of Sp2/0 myeloma cells with spleen cells from BALB/c mice. **Isotype:** IgG₁

Source: Ascites

Immunogen: Native human adiponectin

Format: Purified, Liquid

Purification: >90% pure (SDS-PAGE). Protein A chromatography

Concentration: Lot specific

Affinity Constant: Not determined

Buffer: PBS, pH 7.4

Preservative: 0.1% Sodium azide

Applications: Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Recommended antibody pair for detection of low molecular weight adiponectin in sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
H01290M	H01263M

Storage: Store at 2-8°C.

Warning: This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.